



Selection of an incorrect site is indicated by blinking of the red lamp and buzzing. Successful puncturing at the correct site is indicated by the green lamp.
(Only at Clarke's point and Hochstetter's point)



The epidermis is made of a special material that leaves no needle mark.

Being compact, the device needs a small space for storage.

Cena bez DPH

60.860,00 Kč

Cena s DPH

73.640,60 Kč

Parametry

Kanyly, injekce, napichování - filtr

Hýždě

Množstevní jednotka

ks

- Detailní popis

Napodobení Clarkeova a Hochstetterova bodu. Umožňuje simulaci v pozici na boku.

Kosterní body jako je přední horní kyčelní trn, přední dolní kyčelní trn a trochanter major, mohou být nahmatány na simulátoru pro efektivní nácvik hledání bodů pro zavedení jehly. Může být zavedena tekutina nedporučuje se však používat roztoky léků, ale pro nácvik pouze kohoutkovou vodu.

- **Hmotnost:** Přibližně 1,2 kg
- **Doplňky:** Pásek, odsávací pytlík, úložný obal

Chcete stránku vytisknout?

Skrýt před tiskem nějaké části stránky.



Vytisknout část stránky

Vytisknout celou stránku tak, jak ji vidíte.



Vytisknout celou stránku

(tento box nebude vytištěn)

Úprava stránky před tiskem

Najedte myší na nějaký obsah na stránce, zobrazí se červený obdélník. Kliknutím na něj obsah skryjete. Poté stiskněte tlačítko 'Vytisknout'.



Vytisknout!



Zobrazit vše